

300 evolve[®]_{CT}

Nitrile Powder Free Examination Gloves

3400 Series

- Provides tensile strength and excellent comfort, while not sacrificing tactile sensitivity.
- Ultra 300 Saver Pack reduces storage space and packaging waste.
- Fingertip textured in blue color.
- Meets ASTM D6319 Standards.
- 300 pcs per box (Size XS-L), 250 pcs per box (Size XL) by weight. 10 boxes per carton.

Available Size: XS 3405
 S 3406
 M 3407
 L 3408
 XL 3409



Saver Pack Series Benefits

Maximize Storage Space	Reduce Packaging Waste	Save on Shipping Cost	Reduce Stocking Time
			
			
			

Features	Benefits
Non-Latex	- No latex allergies.
Powder Free	- Eliminates powder allergies and dry skin irritation.
Finger-Textured	- Provides excellent grip.
300pcs Saver Pack	- Triple the quantity of gloves packed in the same box. - Reduces the usage of packaging materials.

Test	Particular	Specification
Dimension (mm) Palm Width	X-Small (3405)	70 ± 10
	Small (3406)	80 ± 10
	Medium (3407)	95 ± 10
	Large (3408)	111 ± 10
	X-Large(3409)	120 ± 10
Dimension (mm) Length	X-Small (3405)	min 220
	Small (3406)	
	Medium (3407)	min 230
	Large (3408)	
	X-Large (3409)	
Thickness (mm)	Palm Thickness	min. 0.05
	Finger Thickness	min. 0.05
Tensile Strength (MPa)	Before Aging	min. 14
	After Aging	min. 14
Elongation (%)	Before Aging	min. 500
	After Aging	min. 400
Powder Residue	Powder Free (mg/glove)	max. 2
Freedom From Holes	G1/AQL 2.5	
Physical Appearance	Rolled beaded cuff Blue color Fingertip-textured surface Ambidextrous	

Remark: The above data represent the norm of the product which may be changed without notice.

Gloves used for protection against chemotherapy drugs exposure should be selected specifically for the type of drugs being used. User is responsible for determining the safe use of glove in chemotherapy treatment. These gloves are for single and transient use only.

Tested for use with Chemotherapy Drugs per ASTM D6978

Tested Chemotherapy Drug & Concentration	Minimum Breakthrough Detection Time
Cisplatin, 1.0 mg/ml	>240 mins
Cyclophosphamide (Cytoxan), 20.0 mg/ml	>240 mins
Dacarbazine (DTIC), 10.0 mg/ml	>240 mins
Doxorubicin Hydrochloride, 2.0 mg/ml	>240 mins
Etoposide (Toposar), 20.0 mg /ml	>240 mins
Fluorouracil, 50.0 mg/ml	>240 mins
Paclitaxel (Taxol), 6.0 mg /ml	>240 mins
Methotrexate, 25 mg/ml	>240 mins
Mitomycin C, 0.5 mg/ml	>240 mins
Vincristine Sulfate, 1.0 mg/ml	>240 mins
Carmustine (BCNU), 3.3 mg/ml	Not recommended
Thiotepa, 10.0 mg/ml	Not recommended

Warning: Do not use with Carmustine and Thiotepa

Packaging Dimension

Inner	Carton
245mmL x 128mmW x 88mmH	460mmL x 266mmW x 265mmH

Intended Application

Cranberry Evolve[®]CT Nitrile Powder Free Examination Gloves is a disposable glove product that is only applicable for examination, medical treatments, and handling contaminated medical devices to prevent cross contamination between users and patients, and not suitable for surgical use.

Handling and Storage

Store in a dry place at room temperature. Exposed product should be shielded from direct sunlight, intense artificial light, x-ray machines and other sources of ozone.

Hazard Identification

If allergic reaction occurs, discontinue use immediately and consult a physician. This product is not made with natural rubber latex.

USGS
BL-03